

H100N Hand Held Pulse Oximeter



The small and lightweight H100N Handheld Pulse Oximeter is designed for continuous monitoring or spot checks in physician offices or clinics, EMS/transport or in acute care settings.

Features:

- LCD display
- RTC (Real Time Clock) display
- Backlight control and automatic power-off function for power saving
- SpO2 and Pulse Rate measurement (Nellcor)
- Temperature measurement
- Display numeric and waveform of SpO2 simultaneously
- Long battery life - up to 50 hours on 4 "AA" or Ni-MH rechargeable batteries
- Powerful data storage capacity
- Data can be transferred to PC for storage, review and printing
- 10 minutes trend graph of SpO2 and Pulse Rate can be reviewed
- Pitch tone
- Audible and visible alarm capability

