

Sonotrax Fetal Doppler Range



Sonotrax II



Sonotrax Basic



Sonotrax Lite

The SonoTrax series of ultrasonic pocket Dopplers are designed to meet the everyday routine requirements for obstetricians, midwives and pregnant women alike.

FEATURES:

- Easy to use
- Ergonomic design
- Waterproof probes
- Compact and light
- Automatic power-off
- Support 2MHZ, 3MHZ, 4MHZ, 8MHZ
- Automatic and manual FHR Counting
- High sensitivity interchangeable probe
- Long time continuous use (approx 8 hrs)
- Accurate FHR detection with clear sound
- Built-in rechargeable poly-lithium battery and AC recharger

Specification

Heart Rate Range:	50—210 bpm
Resolution:	1 bpm
Audio Sampling Frequency:	4KHz / 8KHz
Recording Length:	480 sec / 240 sec
Audio Power Output:	0.5w
Auto Power Off:	1 minute
Size:	32 x 85 x 138mm
Weight:	290g +/- 5g (inc battery)

	FHR Display	Back Light	Probe Auto Detection Display	Battery Status Display	Low Battery Indicator	AC Recharge
Sonotrax Lite	x	x	x	x	P	x
Sonotrax Basic	x	x	P	x	P	x
Sonotrax II	P	P	P	P	P	P



All units in the Sonotrax series have waterproof interchangeable probes.
Optional Extras: Rechargeable Lithium-ion Battery, Charger, Earphones, Carry Bag, 9V Alkaline Battery, AC Adaptor

