

Sonotrax Vascular Ultrasonic Pocket Doppler



The hand held Doppler has been widely adopted as standard equipment for the assessment of peripheral vascular disease. The versatile, inexpensive and portable equipment allows for fast and simple but effective assessment of both arterial venous disease. Besides, it helps physicians to diagnose some diseases such as lower limb and foot ulceration both in diabetics and non diabetics alike.

Sonotrax Vascular offers a traditional 'pencil' style 4 or 8 MHz probe. It can locate both superficial and deep vessels faster and provides stronger signals than other, more expensive equipment.

Features

- Easy to use
- Ergonomic design
- Interchangeable waterproof probes
- Rapid and effective



Specification

Audio Output Power: 0.5W
Auto Shut Off: 1 minute
Size: 32(L) x 85(W) x 138(H) mm
Weight: 290g (including battery)

